

lipomed® bra lipomed® belt

soft self-forming cups - ensuring positive outcomes

The lipomed bra is a multifunctional bra that offers outstanding wearability and ensures that all the specific requirements for the provision of the specific compression following cosmetic breast surgery, are fulfilled.

Properties

- The self-forming cups (cup sizes B + C) do not require additional underwires and prevents the formation of furrows at the lower breast fold.
- Hook-and-eye closure strips on the straps make for quick and easy adjustment to the desired compression and cup size.
- Lateral reinforcements provide additional stability in the lateral portion of the lower breast and also helps shape the breast.
- All closures are finished in such a way that they don't show through under clothing.
- Hypoallergenic, latex-free warp-knit fabric offers an outstanding microclimate and excellent breathability.

lipomed belt

The lipomed belt ensures certain that the implants stay in position after the operation and that your patients achieve the desired result.



Colours

Black

White

Measuring chart lipomed bra

Circumference in cm	S	M	L	XL
Below bust	70 – 80	80 – 90	90 – 105	105 – 125

Measuring chart lipomed belt

Circumference in cm	1	2	3
Below bust	70 – 85	85 – 105	105 – 125

Features

lipomed bra shapes and supports the tissues after cosmetic plastic surgery. Broad straps ensure the highest level of comfort.



lipomed belt

The compression belt prevents dislocation of the implants.



Hook-and-eye closures in the straps

Adjusting compression and cup size is quick and easy.



Combination front closure

The combination of hook-and-eye and zip closures ensures circumferential compression and makes the garment easier to put on.

Importante Information

Patients who notice swelling or sensory disorders whilst wearing the compression garment must consult a doctor immediately. The garments should not be worn for a long period of time unless under medical supervision.

Indications


- Post-operative compression therapy after breast augmentation
- Post-operative compression therapy after breast reduction and lifting
- Compression therapy after mastopexia
- Support for lax connective tissue, e.g. after pregnancy

Contraindications

- Arterial circulatory disturbances
- High risk of thrombosis
- Breathing difficulties
- Intolerability of the material

Important

The garments should not be worn for a long period of time unless under medical supervision.





circaid® Reduction Kit for Trunk and Breast Oedema

Lymphoedema treatment for breast cancer patients

- Adjustable for oedema fluctuation
- Bi-lateral or Uni-lateral mastectomy
- Customisable to patient dimensions

The reduction kit vest is made from inelastic Breathe-O-Prene® material, delivering dynamic compression with superior containment levels. It is trimmed to the patient's size and indication by a trained Certified Lymphoedema Therapist in clinic. Hook and loop enclosures are applied to deliver the necessary compression and secure a comfortable fit.

The reduction kit vest is available in two sizes (small and large) and two lengths to create a perfect fit. The use of the circaid® Reduction Kit leads to better patient independence in the long term.

CRK5L001	Circaid Reduction Kit Vest Long Small
CKR5L003	Circaid Reduction Kit Vest Long Large
CRK5S001	Circaid Reduction Kit Vest Standard Small
CRK5S003	Circaid Reduction Kit Vest Standard Large

Need more information?

✉ sales@medi.au ☎ 1300 04 6334

medi

The vest is designed for use in acute drainage phase Lymphoedema management, to contain the results obtained during treatment and for long-term maintenance. Patients will enjoy the liberating nature of an inelastic garment system they can remove for regular hygiene and re-apply themselves without the aid of a caregiver or clinician.

Product Features

- Trim as needed to fit patient dimensions
- Made from circaid Breath-o-prene® for comfortable containment
- Fully adjustable for oedema fluctuation
- Bi-lateral or uni-lateral mastectomy
- Easy patient application
- Comfortable for day and nighttime wear

Medical Specifications

Indications

- Truncal oedema
- Lymphoedema
- Post-surgical oedema