



THUASNE

MOBIDERM Intimate Bra

Compression and reduction of edema in the breast and chest-wall after breast cancer surgery.

MOBIDERM Intimate Bra

The compressive lymph bra is used for the compression and drainage of edemas in the area of the chest and thorax wall following mammary carcinoma.


The individually adjustable mobiderm pad achieves a local pressure increase on the tissue affected by edemas.



Rear view

Front view

Characteristics

Colors : nude 

Compression value: 7 - 10 mmHg

Product size options

Sizes

Underbust size INT	Size	68-72 70	73-77 75	78-82 80	83-87 85	88-92 90	93-97 95	98-102 100	103-107 105	108-112 110	113-117 115
		XS		S		M		L		XL	
Cup size	AA / A	80 - 89		90 - 99		100 - 109		112 - 119		120 - 129	
	B / C	84 - 93		94 - 103		104 - 113		116 - 123		124 - 133	
	D / E	88 - 97		98 - 107		108 - 117		120 - 127		128 - 137	

To find your correct size

1. **Measure under the bust in CM** and find the corresponding measurement on the above chart *UNDERBUST SIZE INT*.
2. **Measure the fullest part of your breasts** and find the measurement in *CUP SIZE*.

E.G. 1. Measure under bust is 95 cm.

E.G. 2. Cup Size is 110 cm, therefore I would choose MEDIUM B/C.

Discover it
right now!



Or visit
mobiderm.com.au

Indications: Acute or chronic edemas of the breast and/or the breast wall. Postoperative edemas or edemas after radio therapy (only on intact skin). Chronic infiltration after reconstruction, remaining chronic breast wall edemas. The availability of this product may vary from country to country or region to region, due to specific local regulatory registration procedures or conditions of sale. The device referred to in this document is CE marked in accordance with European Regulation 2017/745 concerning medical devices. Carefully read the instructions for use, indications and contraindications of the products. We reserve the right to make changes to our products without notice to the descriptions and illustrations in this document.



Thuasne Deutschland GmbH
Im Steinkamp 12 · D-30938 Burgwedel - Germany
Contact: sales.export@thuasne.fr

Thuasne SAS - SIREN/RCS Nanterre 542 091 186 - capital 1 950 000 euros
120, rue Marius AUFAN 92300 Levallois-Perret (France)