

Comfiwave

Breast Band

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient's skin that help soften the tissues while promoting lymphatic drainage.


Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and even overnight.



Compression

mmHg	Well suited for
15-21	Breast discomfort, mild to severe lymphoedema with or without pitting tissues and shape distortion. Post breast surgery.

Sizing



Circumferences cm	S	M	L	XL
Above breast	70 - 88	80 - 102	90 - 116	100 - 130
Over breast	70 - 92	80 - 108	90 - 126	100 - 144
Below breast	70 - 80	80 - 92	90 - 104	100 - 116

Comfort

Gentle bust coverage

A softer knit around the chest for a comfortably reduced compression level.

Support

Underarm support

Increased elastic thread count under the arms for stronger compression where you need it most.

Use post-operatively

Following breast augmentation or for those who have had breast conserving surgery such as lumpectomy.

“WOW!! This has significantly increased the effective management of my breast lymphoedema. I only wear it overnight and, every morning, my ‘lymph boob’ feels very, very close to just ‘normal’.”

Lymphoedema patient

Codes

	S		M		L		XL	
REGULAR	CW1-TN2		CW1-TN3		CW1-TN4		CW1-TN5	
Unit of Sale	Single							